

Product Datasheet

Human OX-2 protein (orb555807)

Description	Recombinant Human OX-2/MOX1/CD200 (C-Fc)
Endotoxins	Less than 0.1 EU/ug
Preservatives	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-hFc
Note	For research use only
Protein Sequence	Gln31-Gly232
Source	Human 293 Cells
Uniprot ID	P41217
Expiration Date	6 months from date of receipt.