

Product Datasheet

Human ERp44 protein (orb555817)

Description	Recombinant human ER Resident Protein 44/ERp44/TXNDC4 (C-6His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris HCl 10%glycerol,PH7.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris HCl, 10% Glycerol, pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	Glu30-Asp402
Source	Human Cells
Uniprot ID	Q9BS26
Expiration Date	6 months from date of receipt.