

Product Datasheet

Human TCbIR protein (orb555826)

Description	Recombinant human Transcobalamin II Receptor/TCbIR/8D6A/CD320 (C-Fc)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc
Note	For research use only
Protein Sequence	Ser36-Val231
Source	Human Cells
Uniprot ID	Q9NPF0
Expiration Date	6 months from date of receipt.