

Product Datasheet

Human SEMA3A protein (orb555829)

Description	Recombinant human Semaphorin-3A/SEMA3A (C-6His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris 150mM NaCl 1mM EDTA 20%Glycerol, pH 8.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl 1mM EDTA 20% Glycerol, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	Asn21-Val771(Arg552Ala, Arg555Ala)
Source	Human Cells
Uniprot ID	Q14563
Expiration Date	6 months from date of receipt.