

Product Datasheet

Human PLAT protein (orb555831)

Description	Recombinant human Tissue-Type Plasminogen Activator/PLAT (C-6His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM MES,150mM NaCl,0.2mM GaCl2,pH5.5.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mMMES,150mMNaCl,5%trehalose,0.2mMCaCl2,pH5.5.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Protein Sequence	Ser36-Pro562
Source	Human Cells
Uniprot ID	P00750
Expiration Date	6 months from date of receipt.