

Product Datasheet

Mouse DR6 protein (orb555843)

Description	Recombinant mouse Death Receptor 6/DR6/TNFRSF21/CD358 (C-Fc-6His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc&His
Note	For research use only
Protein Sequence	Gln42-His349
Source	Human Cells
Uniprot ID	Q9EPU5
Expiration Date	6 months from date of receipt.