

Product Datasheet

Human LRG1 protein (orb555844)

Description	Recombinant human Leucine-Rich alpha-2-Glycoprotein/LRG1 (C-6His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 20mM NaCl, pH 7.5.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Tris-HCl, 20mM NaCl, pH 7.5.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Protein Sequence	Val36-Gln347
Source	Human Cells
Uniprot ID	AAH34389.1
Expiration Date	6 months from date of receipt.