

Product Datasheet

Human IgG1 Fc protein (orb555846)

Description	Recombinant human IgG1 Fc (Asp239Glu, Leu241Met)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Protein Sequence	Asp104-Lys330(Asp239Glu,Leu241Met)
Source	Human Cells
Uniprot ID	P01857
Expiration Date	6 months from date of receipt.