

Product Datasheet

Human NGAL protein (orb555855)

Description	Recombinant human NGAL/Lipocalin-2/LCN2 (C-6His, human Cells)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of PBS, 50% glycerol,pH7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, 50% Glycerol, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	Gln21-Gly198
Source	Human Cells
Uniprot ID	P80188
Expiration Date	6 months from date of receipt.