

Product Datasheet

SCRIB antibody (orb555861)

Description

Goat polyclonal antibody to SCRIB

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IF, WB

Applications
Target

SCRIB / Scribble

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

WB: Approx 200-250kDa band observed in lysates of cell lines HeLa, HepG2 and U2OS (calculated MW of 178kDa according to NP_874365.3). This molecular weight is routinely observed by other sources. Recommended concentration: 1-3µg/ml. Primary incubation 1 hour . Notes: Immunofluorescence: Strong expression of the protein seen in the cell junctions of MCF7 cells. Recommended concentration: 10µg/ml.

Protein Sequence

PEGPGKEKELPGQ

Clonality

Polyclonal

MW

178; 175

Entrez

23513

Dilution Range

P-ELISA: 1:1000., IF: 10 µg/ml

Expiration Date

12 months from date of receipt.