



Product Datasheet

Rabbit IgG isotype control (PE) [Out of stock] [Out of stock] (orb555865)

Description	Rabbit IgG isotype control (PE)
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% glycerol.
Form/Appearance	Liquid
Concentration	1.0 mg/ml
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
SH-SY5Y
usin...



Flow
cytometric
analysis of
SH-SY5Y
usin...