


Product Datasheet


glyceraldehyde-3-phosphate dehydrogenase Antibody (orb555879)

Description

Rabbit polyclonal antibody to GAPDH


Species/Host

Rabbit


Immunofluorescence
analysis of HeLa
cell...

Reactivity

Bacteria, Bovine, Canine, Drosophila, E. coli, Fish, Gallus, Hamster, Human, Insect, Mammal, Monkey, Mouse, Other, Plant, Porcine, Rabbit, Rat, Yeast, Zebrafish

Conjugation


Unconjugated

Tested

ICC, IHC-P, IP, WB

Applications
Immunogen

Recombinant protein encompassing a sequence within the center region of human GAPDH. The exact sequence is proprietary.


Immunohistochemical
staining of human
sa...

Target


GAPDH

Preservatives

0.025% ProClin 300. 1XPBS, 1% rAlbumin, 20% Glycerol (pH7)

Form/Appearance

Liquid: 1XPBS, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative.


Immunohistochemical
staining of mouse
in...

Concentration

0.69 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: PC-12 , Rat2 , Neuro 2A , C8D30 , NIH-3T3 , Raw264.7 , C2C12 , E.coli M15

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

36

Uniprot ID

P04406

Dilution Range

WB: 1:5000-1:100000 , IF: 1:100-1:1000 , IHC-P: 1:100-1:1000 , IP: 1:100-1:500