

Product Datasheet

GAD65 antibody (orb555899)

Description

Rabbit polyclonal antibody to GAD65

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IHC-Fr, IHC-P, WB

Applications
Immunogen

Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human GAD65. The exact sequence is proprietary.

Target

GAD65

Preservatives

0.01% Thimerosal. 1XPBS (pH7), 20% Glycerol

Form/Appearance

Liquid: 1XPBS, 20% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: mouse brain , rat brain

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW



65

Uniprot ID
[Q05329](https://www.uniprot.org/entry/Q05329)
Dilution Range


WB: 1:1000-1:10000 , IF: 1:100-1:1000 , IHC-P: 1:100-1:1000 , IHC-Fr: 1:100-1:1000

Expiration Date


12 months from date of receipt.



Immunofluorescence analysis of cortical ...



Immunohistochemical staining of RT2 xeno...



Immunocytochemistry analysis of mouse re...