

Product Datasheet

ADORA2A antibody (orb555909)

Description

Rabbit polyclonal antibody to ADORA2A

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human Adenosine A2A-R. The exact sequence is proprietary.

Target

ADORA2A

Preservatives

0.01% Thimerosal. 1XPBS (pH7), 10% Glycerol

Form/Appearance

Liquid: 1XPBS, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

0.5 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: 293T , rat liver

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

45

Uniprot ID


[P29274](#)

Dilution Range


WB: 1:500-1:3000

Expiration Date


12 months from date of receipt.



Western blot analysis of Rat tissue using...



Western blot analysis of 293T cells using...



Western blot analysis of HeLa(Lane 1), M...