

Product Datasheet

galanin receptor 2 Antibody (orb555914)

Description

Rabbit polyclonal antibody to GALR2


Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

ICC, IHC-P, WB

Applications

Immunofluorescence
analysis of
hippocamp...

Immunogen

Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human GALR2. The exact sequence is proprietary.

Target


GALR2

Preservatives

0.01% Thimerosal. 1XPBS (pH7), 50% Glycerol

Form/Appearance

Liquid: 1XPBS, 50% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.


Immunohistochemical
staining of rat fore...

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: A431

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW


42

Uniprot ID
Q43603
Dilution Range

WB: 1:500-1:3000 , IF: 1:100-1:1000 , IHC-P: 1:100-1:1000

Expiration Date

12 months from date of receipt.


Immunohistochemical
staining of Cal27 Xe...