

Product Datasheet

Histamine H4 Receptor antibody (orb555916)

Description

Rabbit polyclonal antibody to Histamine H4 Receptor

kDa 130 A B C

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

Carrier-protein conjugated synthetic peptide encompassing a sequence within the N-terminus region of human Histamine H4 Receptor. The exact sequence is proprietary.

Target

Histamine H4 Receptor

Preservatives

0.01% Thimerosal. 1XPBS (pH7), 50% Glycerol

Form/Appearance

Liquid: 1XPBS, 50% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: Jurkat , Raji , NCI-H929 , human shHRH4-transfected HeLa cells

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

44

Uniprot ID


[Q9H3N8](#)

Dilution Range

WB: 1:500-1:3000

Expiration Date

12 months from date of receipt.



Western blot analysis of Raji cells(Lane...