


Product Datasheet

nuclear receptor subfamily 3 group C member 1 Antibody (orb555939)

Description

Rabbit polyclonal antibody to


Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IHC-P, IP, WB

Applications
Immunogen

Recombinant protein encompassing a sequence within the N-terminus region of human Glucocorticoid Receptor. The exact sequence is proprietary.

Target

Glucocorticoid Receptor

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: A431, mouse heart, rat heart, HeLa

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW


86

Uniprot ID
P04150
Dilution Range


WB: 1:500-1:3000, IHC-P: 1:100-1:1000, IP: 1:500-1:1000

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of HeLa xen...



Immunoprecipitation analysis of HeLa Who...



Western blot analysis of Rat tissue usin...