

Product Datasheet

APP antibody (orb555949)

Description

Rabbit polyclonal antibody to APP

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IHC-Fr, IHC-P, WB

Applications

Immunogen

Recombinant protein encompassing a sequence within the C-terminus region of human APP. The exact sequence is proprietary.

Target

APP

Preservatives

0.01% Thimerosal. 1XPBS (pH7), 10% Glycerol

Form/Appearance

Liquid: 1XPBS, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

0.2 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: HepG2 , U87-MG , SK-N-SH , IMR32 , SK-N-AS

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

87

Uniprot ID



[P05067](#)

Dilution Range


WB: 1:500-1:3000 , IF: 1:100-1:1000 , IHC-P: 1:100-1:1000

Expiration Date


12 months from date of receipt.



Immunofluorescence analysis of SK-N-SH C...



Immunohistochemical staining of human to...



Immunohistochemical staining of Rat brai...