

Product Datasheet

NHERF2 antibody (orb555955)

Description

Rabbit polyclonal antibody to NHERF2

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

Recombinant protein encompassing a sequence within the center region of human NHERF-2. The exact sequence is proprietary.

Target

NHERF2

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

5.91 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: Mouse liver , NT2D1

Isotype

IgG

Clonality

Polyclonal

Purity

Affinity purified by Protein A.

MW

37

Uniprot ID


[Q15599](#)

Dilution Range


WB: 1:500-1:3000

Expiration Date

12 months from date of receipt.



Western blot analysis of Whole cell extr...



Western blot analysis of mouse liver usi...