

Product Datasheet

43715 antibody (orb555964)

Description

Rabbit polyclonal antibody to 43715

Species/Host

Rabbit

Reactivity

Human, Zebrafish

Conjugation

Unconjugated

Tested

IHC, IHC-P, WB

Applications

Immunogen

Recombinant protein encompassing a sequence within the center region of human SEPT7. The exact sequence is proprietary.

Target

43715

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1.02 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: HeLaS3

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

51

Uniprot ID


[Q16181](#)

Dilution Range


WB: 1:500-1:3000 , IHC-P: 1:100-1:1000

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of human ga...



Western blot analysis of HeLa S3 cells u...



Western blot analysis of MOLT4 cells usi...