

Product Datasheet

**acetylcholinesterase (Cartwright blood group)
Antibody (orb555971)**

Description

Rabbit polyclonal antibody to AChE


Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IHC-Fr, IHC-P, WB

Applications

Immunohistochemical staining of rat fore...

Immunogen

Recombinant protein encompassing a sequence within the C-terminus region of human AChE. The exact sequence is proprietary.

Target


AChE

Preservatives

0.025% ProClin 300. 1XPBS pH 7, 20% Glycerol

Form/Appearance

Liquid: 1XPBS, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative.




Immunohistochemical staining of mouse CO...

Concentration

0.73 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.



Immunocytochemistry analysis of Rat brai...

Note

For research use only

Application notes

Positive control: Raji , mouse brain , rat brain

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

68

Uniprot ID

P22303

Dilution Range

WB: 1:500-1:3000 , IF: 1:100-1:1000 , IHC-P: 1:100-1:1000 , IHC-Fr: 1:100-1:1000

Expiration Date

12 months from date of receipt.