

## Product Datasheet

**eukaryotic translation initiation factor 4 gamma 2  
Antibody (orb555977)**

**Description**

Rabbit polyclonal antibody to DAP5


**Species/Host**

Rabbit

**Reactivity**

Human, Mouse

**Conjugation**

Unconjugated

**Tested**

ICC, IHC-P, WB

**Applications**

Immunofluorescence  
analysis of HeLa  
cell...

**Immunogen**


Recombinant protein encompassing a  
sequence within the N-terminus region of  
human DAP5. The exact sequence is  
proprietary.

**Target**

DAP5

**Preservatives**

0.01% Thimerosal. 1XPBS (pH7), 40%  
Glycerol


**Form/Appearance**

Liquid: 1XPBS, 40% Glycerol (pH7). 0.01%  
Thimerosal was added as a preservative.


Immunohistochemical  
staining of mouse  
hi...

**Concentration**

1 mg/ml

**Storage**

Store as concentrated solution.  
Centrifuge briefly prior to opening vial.  
For short-term storage (1-2 weeks), store  
at 4°C. For long-term storage, aliquot and  
store at -20°C or below. Avoid multiple  
freeze-thaw cycles.


**Note**

For research use only

**Application notes**

Positive control: A431 , Mouse brain

Immunohistochemical  
staining of human  
sk...

**Isotype**

IgG

**Clonality**

Polyclonal

**Purity**

Purified by antigen-affinity  
chromatography.

**MW**

102

**Uniprot ID**

**P78344**

**Dilution Range**

WB: 1:500-1:3000 , IF: 1:100-1:1000 ,  
IHC-P: 1:100-1:1000

**Expiration Date**

12 months from date of receipt.