

Product Datasheet

Sec61 translocon gamma subunit Antibody (orb555984)

Description

Rabbit polyclonal antibody to Sec61 gamma

Species/Host

Rabbit

Reactivity

Human, Zebrafish

Conjugation

Unconjugated

Tested

Applications

IHC, WB

Immunogen

Carrier-protein conjugated synthetic peptide encompassing a sequence within the center region of human Sec61 gamma subunit. The exact sequence is proprietary.

Target

Sec61 gamma subunit

Preservatives

0.01% Thimerosal. 1XPBS (pH7), 20% Glycerol

Form/Appearance

Liquid: 1XPBS, 20% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: 293T

Isotype

IgG

Clonality

Polyclonal

Purity

Affinity purified by Protein A.

MW

8

Uniprot ID


P60059

Dilution Range

WB: 1:500-1:3000

Expiration Date

12 months from date of receipt.



Western blot analysis of 293T cells using...