


Product Datasheet


**DnaJ heat shock protein family (Hsp40) member
B6 Antibody (orb556002)**

Description

Rabbit polyclonal antibody to DNAJB6



Immunofluorescence
analysis of HeLa
cell...



Western blot
analysis of Whole
cell extr...

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Recombinant protein encompassing a sequence within the center region of human DNAJB6. The exact sequence is proprietary.
Target	DNAJB6
Preservatives	0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol
Form/Appearance	Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Concentration	0.46 mg/ml
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Note	For research use only
Application notes	Positive control: MDA-MB-231
Isotype	IgG
Clonality	Polyclonal
Purity	Purified by antigen-affinity chromatography.
MW	36
Uniprot ID	O75190
Dilution Range	WB: 1:500-1:3000 , IF: 1:100-1:1000
Expiration Date	12 months from date of receipt.