

Product Datasheet

ADAMTSL2 antibody (orb556015)

Description

Rabbit polyclonal antibody to ADAMTSL2

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

Recombinant protein encompassing a sequence within the C-terminus region of human ADAMTSL2. The exact sequence is proprietary.

Target

ADAMTSL2

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: H1299

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

105

Uniprot ID


[Q86TH1](#)

Dilution Range

WB: 1:500-1:3000 , IHC-P: 1:100-1:1000


Expiration Date

12 months from date of receipt.



Immunohistochemical staining of SkHep1

X...



Western blot analysis of H1299 cells using...