



Product Datasheet

NDUFV1 antibody (orb556026)


Description	Rabbit polyclonal antibody to NDUFV1
Species/Host	Rabbit
Reactivity	Human, Mouse, Zebrafish
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant protein encompassing a sequence within the center region of human NDUFV1. The exact sequence is proprietary.
Target	NDUFV1
Preservatives	0.025% ProClin 300. 1XPBS (pH7), 20% Glycerol
Form/Appearance	Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Concentration	1 mg/ml
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Purified by antigen-affinity chromatography.
MW	51
Uniprot ID	P49821
Dilution Range	WB: 1:500-1:3000 , IHC-P: 1:100-1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse kidney tissue.



Immunohistochemical staining of mouse liver tissue.



Immunohistochemical staining of U373 xenograft tissue.