

Product Datasheet


alpha Glucosidase II antibody (orb556031)

Description


Rabbit polyclonal antibody to alpha Glucosidase

GTX102237


Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	Recombinant protein encompassing a sequence within the N-terminus region of human alpha Glucosidase II. The exact sequence is proprietary.
Target	alpha Glucosidase II
Preservatives	0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol
Form/Appearance	Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Concentration	1 mg/ml
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Note	For research use only
Application notes	Positive control: Mouse brain , 293T , A431 , H1299 , HeLaS3 , HepG2 , Molt-4 , Raji , HeLa
Isotype	IgG
Clonality	Polyclonal
Purity	Purified by antigen-affinity chromatography.
MW	107
Uniprot ID	Q14697
Dilution Range	WB: 1:500-1:3000 , IF: 1:100-1:1000 , IHC-P: 1:100-1:1000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of SW480 Xe...



Immunohistochemical staining of human br...