

Product Datasheet

NADH:ubiquinone oxidoreductase core subunit V2 Antibody (orb556032)

Description

Rabbit polyclonal antibody to NDUFV2



Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunohistochemical staining of mouse heart tissue

Immunogen

Recombinant protein encompassing a sequence within the center region of human NDUFV2. The exact sequence is proprietary.

Target


NDUFV2

Preservatives

0.01% Thimerosal. 1XPBS (pH7), 10% Glycerol

Form/Appearance

Liquid: 1XPBS, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.



Immunohistochemical staining of mouse heart tissue

Concentration

0.44 mg/ml

Storage


Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: mouse liver, rat liver, A431



Western blot analysis of A431 cells using the antibody

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

27

Uniprot ID

[P19404](https://www.uniprot.org/entry/P19404)

Dilution Range

WB: 1:500-1:3000, IHC-P: 1:100-1:1000

Expiration Date

12 months from date of receipt.