

Product Datasheet

KCC4 antibody (orb556050)

Description

Rabbit polyclonal antibody to KCC4

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

Recombinant protein encompassing a sequence within the C-terminus region of human KCC4. The exact sequence is proprietary.

Target

KCC4

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7).
0.01% Thimerosal was added as a preservative.

Concentration

4.71 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: H1299

Isotype

IgG

Clonality

Polyclonal

Purity

Affinity purified by Protein A.

MW


119

Uniprot ID
[Q9Y666](https://www.uniprot.org/entry/Q9Y666)
Dilution Range


WB: 1:500-1:3000

Expiration Date

12 months from date of receipt.



Western blot analysis of H1299 cells using anti-KCC4 antibody.



Western blot analysis of HeLa S3 cells using anti-KCC4 antibody.