


Product Datasheet


DNA polymerase delta 1, catalytic subunit
Antibody (orb556060)

Description


Rabbit polyclonal antibody to DNA



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of human lu...



Western blot analysis of Various whole c...

Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	Recombinant protein encompassing a sequence within the center region of human DNA polymerase delta. The exact sequence is proprietary.
Target	DNA polymerase delta
Preservatives	0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol
Form/Appearance	Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Concentration	0.56 mg/ml
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Note	For research use only
Application notes	Positive control: Rat2 , Molt-4 , 293T , A431 , HeLa , HepG2
Isotype	IgG
Clonality	Polyclonal
Purity	Purified by antigen-affinity chromatography.
MW	124
Uniprot ID	P28340
Dilution Range	WB: 1:500-1:3000 , IF: 1:100-1:1000 , IHC-P: 1:100-1:1000
Expiration Date	12 months from date of receipt.