

Product Datasheet

mitochondrial ribosomal protein S5 Antibody (orb556125)

Description

Rabbit polyclonal antibody to MRPS5



Species/Host

Rabbit

Reactivity

Human, Mouse


Conjugation

Unconjugated

Tested

ICC, IHC-P, WB

Applications



Immunofluorescence analysis of HeLa cell...

Immunogen

Recombinant protein encompassing a sequence within the center region of human MRPS5. The exact sequence is proprietary.

Target


MRPS5

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.



Immunohistochemical staining of mouse ki...

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: 293T, A431, HeLa, HepG2

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

48

Uniprot ID


[P82675](https://www.uniprot.org/entry/P82675)

Dilution Range

WB: 1:500-1:3000, IF: 1:100-1:1000, IHC-P: 1:100-1:1000

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of AGS xeno...