

Product Datasheet

ALAS-E antibody (orb556149)

Description

Rabbit polyclonal antibody to ALAS-E

Species/Host

Rabbit

Reactivity

Human, Mouse, Zebrafish

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

Recombinant protein encompassing a sequence within the center region of human ALAS-E. The exact sequence is proprietary.

Target

ALAS-E

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: K562 , Mouse liver

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

65

Uniprot ID


P22557

Dilution Range


WB: 1:500-1:3000

Expiration Date

12 months from date of receipt.



Western blot analysis of mouse liver using...



Western blot analysis of K562 cell using...