

Product Datasheet

junctional adhesion molecule 2 Antibody (orb556173)

Description

Rabbit polyclonal antibody to JAM-B


Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC, IHC-P, WB

Applications

Immunohistochemical
staining of human
CO...

Immunogen


Recombinant protein encompassing a
sequence within the center region of human
JAM-B. The exact sequence is proprietary.

Target

JAM-B

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7),
10% Glycerol



Immunohistochemical
staining of mouse
sm...

Form/Appearance


Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol
(pH7). 0.01% Thimerosal was added as a
preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge
briefly prior to opening vial. For short-term
storage (1-2 weeks), store at 4°C. For long-
term storage, aliquot and store at -20°C or
below. Avoid multiple freeze-thaw cycles.



Immunohistochemical
staining of TOV-
112D...

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

33

Uniprot ID

P57087

Dilution Range

WB: 1:500-1:3000 , IHC-P: 1:100-1:1000

Expiration Date

12 months from date of receipt.