

Product Datasheet

matrix metalloproteinase 2 Antibody (orb556175)

Description

Rabbit polyclonal antibody to MMP2

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ICC, IHC-Fr, IHC-P, IP, WB

Applications
Immunogen

Recombinant protein encompassing a sequence within the center region of human MMP2. The exact sequence is proprietary.

Target

MMP2

Preservatives

0.025% ProClin 300, 1XPBS pH 7, 1% rAlbumin, 20% Glycerol

Form/Appearance

Liquid: 1XPBS, 1% rAlbumin, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative.

Concentration

0.17 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

74

Uniprot ID

P08253

Dilution Range


WB: 1:500-1:3000, IF: 1:100-1:1000, IHC-P: 1:100-1:1000, IP: 1:100-1:500

Expiration Date


12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of human en...



Immunohistochemical staining of human ov...