


Product Datasheet

SDHB antibody (orb556185)

Description

Rabbit polyclonal antibody to SDHB


Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ICC, IHC-P, WB

Applications

Immunofluorescence
analysis of HeLa
cell...

Immunogen

Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human SDHB. The exact sequence is proprietary.

Target


SDHB

Preservatives

0.01% Thimerosal. 1XPBS (pH7), 40% Glycerol

Form/Appearance

Liquid: 1XPBS, 40% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.


Immunohistochemical
staining of human
CO...

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: THP-1 , NCI-H929 , Neuro2A

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

32

Uniprot ID


P21912

Dilution Range

WB: 1:500-1:3000 , IF: 1:100-1:1000 , IHC-P: 1:100-1:1000

Expiration Date

12 months from date of receipt.


Immunohistochemical
staining of mouse
mu...