

## Product Datasheet

### Alpha B Crystallin ELISA Kit [Out of stock] (orb55621)

|                            |                                                 |
|----------------------------|-------------------------------------------------|
| <b>Description</b>         | Colorimetric detection of<br>alpha B Crystallin |
| <b>Range</b>               | 0.352 - 22.5 ng/ml                              |
| <b>Tested Applications</b> | ELISA                                           |
| <b>Target</b>              | alpha B Crystallin                              |
| <b>Storage</b>             | 4°C                                             |
| <b>Note</b>                | For research use only                           |
| <b>Sample Types</b>        | Cell Lysates   Tissue  <br>Serum                |
| <b>Sensitivity</b>         | 0.009 ng/ml                                     |
| <b>Expiration Date</b>     | Please enquire.                                 |