

## Product Datasheet

### CTD small phosphatase 2 Antibody (orb556218)

**Description**

Rabbit polyclonal antibody to CTDSP2

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested**

WB

**Applications**

**Immunogen**

Recombinant protein encompassing a sequence within the center region of human CTDSP2. The exact sequence is proprietary.

**Target**

CTDSP2

**Preservatives**

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

**Form/Appearance**

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

**Concentration**

0.96 mg/ml

**Storage**

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

**Note**

For research use only

**Application notes**

Positive control: 293T , Mouse brain , PC-12

**Isotype**

IgG

**Clonality**

Polyclonal

**Purity**

Purified by antigen-affinity chromatography.

**MW**

31

**Uniprot ID**


**O14595**

**Dilution Range**


WB: 1:500-1:3000

**Expiration Date**


12 months from date of receipt.



Western blot analysis of PC-12 whole cell...



Western blot analysis of mouse brain usi...



Western blot analysis of 293T cells usin...