

Product Datasheet

PDSS2 antibody (orb556220)

Description

Rabbit polyclonal antibody to PDSS2

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

Recombinant protein encompassing a sequence within the center region of human PDSS2. The exact sequence is proprietary.

Target

PDSS2

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

0.54 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: Mouse brain , A431

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

44

Uniprot ID


Q86YH6

Dilution Range


WB: 1:500-1:3000

Expiration Date

12 months from date of receipt.



Western blot analysis of A431 cells using...



Western blot analysis of mouse brain using...