

Product Datasheet

IPMK antibody (orb556222)

Description

Rabbit polyclonal antibody to IPMK

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ICC, IHC-P, WB

Applications

Immunogen

Recombinant protein encompassing a sequence within the center region of human IPMK. The exact sequence is proprietary.

Target

IPMK

Preservatives

0.025% ProClin 300. 1XPBS pH 7, 20% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

0.69 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: Mouse brain

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

47

Uniprot ID

[Q8NFU5](#)

Dilution Range


WB: 1:500-1:3000 , IF: 1:100-1:1000

Expiration Date


12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Western blot analysis of H1299 cells usi...



Western blot analysis of mouse brain usi...