

Product Datasheet

IL17 Receptor C antibody (orb556252)

Description

Rabbit polyclonal antibody to IL17 Receptor C

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

Recombinant protein encompassing a sequence within the center region of human IL17 receptor C. The exact sequence is proprietary.

Target

IL17 Receptor C

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: HeLaS3

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW


86

Uniprot ID
Q8NAC3
Dilution Range

WB: 1:500-1:3000

Expiration Date

12 months from date of receipt.



Western blot analysis of HeLa S3 cells u...