

## Product Datasheet

### Interferon alpha 8 antibody (orb556267)

## Description

Rabbit polyclonal antibody to Interferon alpha 8

## Species/Host

Rabbit

## Reactivity

Human

## Conjugation

Unconjugated

## Tested

IHC-P, WB

## Applications

## Immunogen

Recombinant protein encompassing a sequence within the center region of human interferon alpha 8. The exact sequence is proprietary.

## Target

Interferon alpha 8

## Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

## Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

## Concentration

1 mg/ml

## Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

## Note

For research use only

## Application notes

Positive control: Raji

## Isotype

IgG

## Clonality

Polyclonal

## Purity

Purified by antigen-affinity chromatography.

## MW

22

## Uniprot ID


[P32881](#)

## Dilution Range


WB: 1:500-1:3000 , IHC-P: 1:100-1:1000

## Expiration Date

12 months from date of receipt.



Immunohistochemical staining of OVCA xen...



Western blot analysis of Raji cells usin...