


Product Datasheet

WBSCR22 antibody (orb556287)

Description

Rabbit polyclonal antibody to WBSCR22


Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ICC, IHC-P, IP, WB

Applications

Immunofluorescence
analysis of HeLa
cell...

Immunogen

Recombinant protein encompassing a sequence within the center region of human WBSCR22. The exact sequence is proprietary.

Target


WBSCR22

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.


Immunohistochemical
staining of mouse
ki...

Concentration

0.42 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: HeLa , HCT116

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW


32

Uniprot ID
O43709
Dilution Range

WB: 1:500-1:3000 , IF: 1:100-1:1000 , IHC-P: 1:100-1:1000 , IP: 1:500-1:1000

Expiration Date

12 months from date of receipt.


Immunohistochemical
staining of Ca922
xe...