

Product Datasheet

malate dehydrogenase 2 Antibody (orb556289)

Description

Rabbit polyclonal antibody to MDH2


Species/Host

Rabbit

Reactivity


Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IHC-P, WB

Applications

Immunofluorescence
analysis of HeLa
cell...

Immunogen

Recombinant protein encompassing a sequence within the center region of human MDH2. The exact sequence is proprietary.

Target


MDH2

Preservatives

0.025% ProClin 300, 1XPBS (pH7), 1% rAlbumin, 20% Glycerol

Form/Appearance

Liquid: 1XPBS, 1% rAlbumin, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative.


Immunohistochemical
staining of mouse
he...

Concentration

0.38 mg/ml

Storage


Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: Mouse brain, rat lung, 293T, A431, HeLa, HepG2


Immunohistochemical
staining of mouse
ki...

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

36

Uniprot ID

P40926

Dilution Range

WB: 1:500-1:20000, IF: 1:100-1:1000, IHC-P: 1:100-1:1000

Expiration Date

12 months from date of receipt.