

## Product Datasheet

### KPNA4 antibody (orb556336)

## Description

Rabbit polyclonal antibody to KPNA4

## Species/Host

Rabbit

## Reactivity

Human, Mouse, Rat

## Conjugation

Unconjugated

## Tested

ICC, IHC-P, WB

## Applications

## Immunogen

Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human KPNA4. The exact sequence is proprietary.

## Target

KPNA4

## Preservatives

0.01% Thimerosal. 1XPBS (pH7), 1% rAlbumin, 20% Glycerol

## Form/Appearance

Liquid: 1XPBS, 1% rAlbumin, 20% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

## Concentration

1.04 mg/ml

## Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

## Note

For research use only

## Application notes

Positive control: 293T, HepG2, Molt-4, Rat2, Neuro2A, GL261, C8D30, NIH-3T3, BCL-1, HeLa

## Isotype

IgG

## Clonality

Polyclonal

## Purity

Purified by antigen-affinity chromatography.

## MW

58

## Uniprot ID


[O00629](https://www.uniprot.org/entry/O00629)

## Dilution Range


WB: 1:500-1:3000, IF: 1:100-1:1000, IHC-P: 1:100-1:1000

## Expiration Date


12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of H1299 xe...



Western blot analysis of 293T(Lane 1), H...