


Product Datasheet

interleukin 1 receptor antagonist Antibody (orb556363)

Description

Rabbit polyclonal antibody to IL1 Receptor



Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IHC-P, WB

Applications

Immunofluorescence
analysis of HeLa
cell...

Immunogen

Recombinant protein encompassing a sequence within the center region of human IL1 Receptor antagonist. The exact sequence is proprietary.

Target

IL1 Receptor antagonist

Preservatives

0.025% ProClin 300, 1XPBS (pH7), 20% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: Neuro2A whole, C8D30, NIH-3T3, Raw264.7, C2C12, PC-12

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

20

Uniprot ID


[P18510](#)

Dilution Range


WB: 1:500-1:3000, IF: 1:100-1:1000, IHC-P: 1:100-1:1000

Expiration Date

12 months from date of receipt.



Immunohistochemical
staining of PC3
xeno...



Western blot analysis
of PC-12 whole cel...