

Product Datasheet

PACAP antibody (orb556384)

Description

Rabbit polyclonal antibody to PACAP

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IHC-Fr, IHC-P, WB

Applications
Immunogen

Recombinant protein encompassing a sequence within the center region of human PACAP. The exact sequence is proprietary.

Target

PACAP

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration

1.02 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: Molt-4 , Raji , mouse spleen

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW

21

Uniprot ID
[Q8WU39](https://www.uniprot.org/entry/Q8WU39)
Dilution Range


WB: 1:500-1:3000 , IHC-P: 1:100-1:1000 , IHC-Fr: 1:100-1:1000

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of human no...



Immunocytochemistry analysis of mouse br...



Western blot analysis of MOLT4(Lane 1), ...