



Product Datasheet

KSR antibody (orb556462)

Description	Rabbit polyclonal antibody to KSR
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein encompassing a sequence within the center region of human Dyrk3. The exact sequence is proprietary.
Target	KSR
Preservatives	0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol
Form/Appearance	Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Concentration	0.68 mg/ml
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Purified by antigen-affinity chromatography.
MW	66
Uniprot ID	O43781
Dilution Range	WB: 1:500-1:3000
Expiration Date	12 months from date of receipt.



Western blot analysis of Raji cells using KSR antibody.



Western blot analysis of mouse liver tissue using KSR antibody.