

## Product Datasheet

### EPM2AIP1 antibody (orb556471)

## Description

Rabbit polyclonal antibody to EPM2AIP1

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

WB

### Applications

### Immunogen

Carrier-protein conjugated synthetic peptide encompassing a sequence within the N-terminus region of human EPM2AIP1. The exact sequence is proprietary.

### Target

EPM2AIP1

### Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

### Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

### Concentration

1 mg/ml

### Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

### Note

For research use only

### Application notes

Positive control: 293T , A431 , H1299

### Isotype

IgG

### Clonality

Polyclonal

### Purity

Purified by antigen-affinity chromatography.

### MW

70

### Uniprot ID


[Q7L775](#)

### Dilution Range


WB: 1:500-1:3000

### Expiration Date

12 months from date of receipt.



Western blot analysis of 293T (Lane 1), A...



Western blot analysis of A549 cells usin...