

Product Datasheet

glutathione S-transferase alpha 1 Antibody (orb556476)

Description

Rabbit polyclonal antibody to GSTA1



Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ICC, IHC-P, WB

Applications



Immunofluorescence analysis of Hep3B cells...

Immunogen

Recombinant protein corresponding to human GSTA1

Target


GSTA1

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.



Immunohistochemical staining of human colon tissue...

Concentration

6 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: mouse liver, HepG2

Isotype

IgG

Clonality

Polyclonal

Purity

Affinity purified by Protein A.

MW

26

Uniprot ID


[P08263](#)

Dilution Range

WB: 1:1000-1:10000, IF: 1:100-1:1000, IHC-P: 1:100-1:1000

Expiration Date

12 months from date of receipt.



Western blot analysis of Hep G2 cells using this antibody.