

Product Datasheet

RAMP1 antibody (orb556540)

Description

Rabbit polyclonal antibody to RAMP1

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human RAMP1. The exact sequence is proprietary.

Target

RAMP1

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid: 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7).
0.01% Thimerosal was added as a preservative.

Concentration

0.87 mg/ml

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Positive control: HepG2 , mouse brain , rat brain

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-affinity chromatography.

MW


17

Uniprot ID
[O60894](https://www.uniprot.org/entry/O60894)
Dilution Range


WB: 1:500-1:3000

Expiration Date

12 months from date of receipt.



Western blot analysis of Various tissue ...



Western blot analysis of Hep G2 cells us...